

**Gas Chromatograph Mass  
Spectrometer  
GCMS-QP2010 Series  
Operation Guide  
For GCMSsolution Ver. 2.6**

**Read the instruction manual thoroughly before you use the product.  
Keep this instruction manual for future reference.**

 **SHIMADZU CORPORATION**  
KYOTO JAPAN

ANALYTICAL & MEASURING INSTRUMENTS DIVISION

This page is intentionally left blank.

In order to ensure the safe use of this product, observe the following points.

- ◆ Follow the procedures described in the instruction manual.
- ◆ Observe precautionary information.
- ◆ Do not disassemble or modify this product without permission.
- ◆ If repairs are required, contact your Shimadzu representative.

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries. Other product and company names mentioned herein are trademarks or registered trademarks of their respective companies.

"Microsoft® Windows®" is abbreviated to "Windows".

<< Notices >>

- All rights are reserved, including those to reproduce this document or parts thereof in any form without permission in writing from Shimadzu Corporation.
- Information in this document is subject to change without notice and does not represent a commitment on the part of the vendor.
- Any errors or omissions which may have occurred in this document despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection. We appreciate notification of any errors or omissions.
- Shimadzu Corporation is not responsible for errors or injuries resulting from following the instructions in this document.
- Shimadzu Corporation is not responsible for errors or injuries resulting from the use of equipment by the customer.
- Replacement parts for this product will be available for a period of seven (7) years after discontinuation of the product. Thereafter, such parts may cease to be available. Note, however, that the availability of parts not manufactured by Shimadzu Corporation is determined by the manufacturer of those parts.
- The contents of PC hard drives may be lost due to unforeseen circumstances. In order to protect important data, be sure to make backup copies on a regular basis.

© 2010 Shimadzu Corporation. All rights reserved.



## Introduction

Thank you for purchasing the GCMS-QP2010-series gas chromatograph mass spectrometer. This document is intended to explain basic operations to first-time users. Refer to the instruction manual or the appendix of this document for more details, including information related to maintenance and inspection.

This document assumes that the user has a working knowledge of Windows. There are references to functions and terminology specific to Windows; refer to a Windows user manual as necessary. Users who have never worked with Windows should read a Windows user manual before using this document.

# Product Warranty

Shimadzu Corporation provides the following warranty for this product.

## Details

1. Period: Please contact your Shimadzu representative for information concerning the extent of your product's warranty.
2. Terms: If a fault attributable to Shimadzu Corporation occurs within the warranty period, we will perform the necessary repairs or part replacement free of charge. Note, however, that it may not be possible to replace products that have short life cycles, such as PCs and their peripheral devices and parts, with the same models.
3. Limitation of Liability
  - 1) In no event will Shimadzu be liable for any lost revenue, profit or data, or for special, indirect, consequential, incidental or punitive damages, however caused regardless of the theory of liability, arising out of or related to the use of or inability to use the product, even if Shimadzu has been advised of the possibility of such damage.
  - 2) In no event will Shimadzu's liability to you, whether in contract, tort (including negligence), or otherwise, exceed the amount you paid for the product.
4. Items Not Covered: The following items are not covered by this warranty, even if they occur during the warranty period.
  - 1) Faults resulting from incorrect use
  - 2) Faults resulting from repairs or modifications implemented by parties other than Shimadzu Corporation or companies designated by Shimadzu Corporation
  - 3) Faults occurring after use in combination with hardware or software other than that designated by Shimadzu Corporation
  - 4) Faults in equipment and damage to data and software, including the OS, resulting from computer viruses
  - 5) Faults in equipment and damage to data and software, including the OS, resulting from power failures, including power interruptions and momentary voltage drops
  - 6) Faults in equipment and damage to data and software, including the OS, resulting from turning OFF the power switch on the equipment without following the proper shutdown procedures
  - 7) Faults not originating in the equipment itself
  - 8) Faults occurring after use in severe environmental conditions, such as those subject to high temperatures, high humidity levels, corrosive gases, or vibrations
  - 9) Faults resulting from fires, earthquakes and other natural disasters, contamination by radioactive or toxic substances, wars, riots, criminal activities, and other types of force majeure
  - 10) Faults occurring after the product is moved or transported following initial installation
  - 11) Faults occurring in consumable parts or parts dependent on them





Note: Recording media, such as floppy disks and CD-ROMs, are also regarded as consumable parts.

\* If a warranty certificate or a similar form of document is provided with the product, or if a contract specifying warranty terms has been signed, then the regulations specified in such a document take precedence.  
The warranty period for products with special specifications or for system products is specified separately.

## About This Operation Guide

### Notation

This operation guide uses the notation described below.

Notation	Meaning
 <b>WARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in serious injury or possibly death.
 <b>CAUTION</b>	Indicates a potentially hazardous situation which, if not avoided, may result in minor to moderate injury or equipment damage.
 <b>NOTE</b>	Indicates additional information that is provided to ensure the proper use of this product.
 <b>Reference</b>	Indicates the location of related information. "Items indicated by this notation mean that this (new) feature/function is downward compatible with version 2.5. The new features are designed to make the software easier to use and/or increase productivity. However, the user still has the option to either use the same method used in version 2.5 or use the new features by following the steps provided in this manual."
[ ]	Indicates items displayed on the screen, such as buttons, menu selections, settings, windows, and icons. Example: Click [OK].

## Safety Precautions

These safety precautions contain important safety-related information. Be sure to read them before using the product and observe them during use.

### WARNING

- Making internal repairs to the product is potentially dangerous. Such repairs must be performed by specially trained Shimadzu personnel.
- Do not disassemble or modify the product without authorization.  
Doing so may compromise safety.
- Read the instruction manual thoroughly before handling or operating the equipment, and be sure to following the procedures described.  
Not handling the equipment as described is potentially dangerous.

### ■ Installation Site Precautions

#### WARNING

- The solvents used with the gas chromatograph mass spectrometer may be flammable or toxic. Install the product in a well-ventilated room.  
Otherwise, solvent vapors may cause poisoning, or ignite and cause a fire.
- Do not use this product in an environment containing combustible gases.  
Doing so may cause a fire.
- Do not place flammable materials near the column oven exhaust at the back of the gas chromatograph.  
Doing so may cause a fire.
- Install the product on a surface that is level, stable, and strong enough to support the product's weight.  
Otherwise, the product may tip over or fall off the surface.



#### NOTE

Do not install the product in locations subject to large amounts of corrosive gas or dust. Otherwise, performance may be adversely affected and the product's service life may be shortened.

## ■ High-Pressure Gas Precautions

### WARNING

- High-pressure gas cylinders are used to supply the carrier gas. Follow the instructions received from the cylinder suppliers and handle the cylinders carefully.

Not doing so may cause poisoning or a fire.

- Keep the cylinders in a well-ventilated outdoor location that is not directly exposed to sunlight, and use pipes to convey the gas indoors. For liquefied gases, this is required by law.
- Ensure that the temperature of gas cylinders never exceeds 40 °C and that there are no naked flames within 2 m of the cylinders.
- Ensure that the installation site is well ventilated and, as part of the daily inspection procedure, check for gas leaks using soapy water. Do not smoke or use open flames within 5 m of equipment using highly combustible gases, such as acetylene, hydrogen, and propane, or potentially combustible gases, such as oxygen and nitrous oxide. Keep a suitable fire extinguisher nearby at all times.
- Secure the cylinders with ropes or by some other method to prevent them from falling over.
- Be sure to use oil-free pressure reducing valves. Also, do not use pipes in which oil is present on the inner surfaces that make contact with the gas.
- After using the gas, close the main valve immediately.
- Check that the pressure gauges are functional at least once every three months.
- Warning signs (adhesive aluminum plates) that indicate hydrogen gas use are provided free of charge. Contact your Shimadzu representative in cases of particular necessity.

Legal authorization is required to use cylinders with a capacity of 300 m<sup>3</sup> or greater.

Refer to high-pressure gas control laws, liquid petroleum gas safety regulations, general high-pressure gas safety regulations, and fire safety laws for more information.

## ■ Operation Precautions

### WARNING

- Wear safety glasses when handling solvents or when injecting samples into the gas chromatograph. If solvent gets into the eyes, it may cause blindness. If solvent does get into the eyes, immediately flush with large amounts of water and seek medical attention.
- Do not place solvents near PCs, printers, or other types of office equipment. Doing so may cause a fire or equipment failure.
- Do not use flammable sprays (e.g., hair sprays and insecticide sprays) near the product. They may ignite and cause a fire.



## Handling Emergencies

Take the measures described below in the event of an emergency, such as a malfunction of the gas chromatograph mass spectrometer.

Before resuming operation, take appropriate precautions and, if necessary, contact your Shimadzu representative.

### ■ Emergency Shutdown Procedure

- 1 Turn OFF the gas chromatograph mass spectrometer.
- 2 Turn OFF all accessories.
- 3 Close the main valves for the pipes supplying carrier gas, hydrogen, and air.
- 4 Disconnect the power supply.
  - ◆ If the power cable is attached to a switchboard, turn OFF the switchboard.
  - ◆ If the power cable is plugged into an outlet, unplug the cable.



This page is intentionally left blank.



# Contents

---

## **1 Overall Configuration of GCMSsolution**

---

1.1	Programs .....	1
1.2	Flowchart of Operating Procedure .....	2

## **2 Starting GC/MS**

---

2.1	Turning ON the Power .....	3
2.2	Layout of Operating Areas .....	4
2.3	Inspecting Consumable Items and Maintenance Parts .....	5
2.4	System Configuration .....	6
2.4.1	Setting the Modules Used for Analysis .....	6
2.4.2	Checking Column Information .....	7
2.4.3	Enabling the Modules Used for Analysis .....	9
2.5	Vacuum System Startup .....	10
2.6	Checking for Vacuum Leakage .....	11
2.7	Autotuning .....	13
2.7.1	Setting Analysis Conditions .....	13
2.7.2	Executing Autotuning .....	14
2.7.3	Checking Autotuning Results .....	17

## **3 Creating a Folder**




---

3.1	Creating a Folder with Data Explorer .....	18
-----	--	----

## **4 Qualitative Analysis**


---

4.1	Creating a Method File .....	20
4.1.1	Setting Autosampler Parameters .....	21
4.1.2	Setting GC Parameters .....	21
4.1.3	Setting MS Parameters .....	22
4.1.4	Setting Similarity Search Parameters .....	23
4.1.5	Saving the Method File .....	25

<b>4.2</b>	Repeating Autotuning .....	25
<b>4.3</b>	Sequential Analysis .....	26
4.3.1	Creating a Batch File.....	26
4.3.2	Editing a Batch File  Ver. 2.5 .....	29
4.3.3	Saving Batch Files .....	32
4.3.4	Executing Sequential Analysis .....	33
<b>4.4</b>	Analyzing Data .....	34
4.4.1	Loading Data Files .....	34
4.4.2	Zoom in on a Peak in the Chromatogram and Displaying Mass Spectra .....	35
4.4.3	Removing the Background.....	37
4.4.4	Searching for Similarity  Ver. 2.5 .....	38
4.4.5	Displaying Mass Chromatograms .....	39
4.4.6	Registering Spectra Displayed for Target Compounds  Ver. 2.5 .....	40
4.4.7	Editing the Spectrum Process Table .....	40
4.4.8	Saving Data Files .....	42
<b>4.5</b>	Printing Qualitative Analysis Reports .....	43
4.5.1	Loading Report Formats .....	43
4.5.2	Editing Report Formats .....	44
4.5.3	Outputting Reports .....	47

## 5 Quantitative Analysis

---

<b>5.1</b>	Creating a Method File .....	48
5.1.1	Creating a Compound Table .....	48
5.1.2	Creating a SIM Table .....	53
<b>5.2</b>	Sequential Analysis .....	55
5.2.1	Creating a Batch File.....	55
5.2.2	Editing a Batch File  Ver. 2.5.....	58
5.2.3	Saving Batch Files .....	58
5.2.4	Executing Sequential Analysis .....	59
<b>5.3</b>	Analyzing Data .....	60
5.3.1	Checking and Correcting Calibration Curves .....	60
5.3.2	Re-quantifying after Correcting a Calibration Curve .....	67
5.3.3	Checking and Correcting Quantitation Results .....	69
<b>5.4</b>	Printing Quantitative Analysis Reports .....	73
5.4.1	Creating and Outputting Quantitative Analysis Reports .....	73

## 6 Shutting Down GC/MS

---

- 6.1 Vacuum System Shutdown.....76
- 6.2 Turning OFF the Power .....77

## Appendix A File Format

---

## Appendix B Viewing Help

---


- B.1 Viewing Help .....79
  - B.1.1 Displaying Help from the Assistant Bar .....79
  - B.1.2 Displaying Help from the Menu Bar .....79
  - B.1.3 Displaying Help with the F1 Key .....80

## Appendix C Using Data Explorer

---

## Appendix D Maintenance

---

- D.1 Maintenance .....83
- D.2 Easy sTop  
(Applicable to QP2010 Ultra and QP2010 SE models) .....85
- D.3 Reset Procedure for Usage Frequencies and Usage Times  Ver.2.5 .....87
- D.4 Changing Replacement Guidelines for septa and Glass Inserts .....89

## Appendix E Single Analysis (Manual Injection)

---

## Appendix F Index Searches

---

## Appendix G Displaying Chromatograms

---

- G.1 Displaying Chromatograms from Fragment Tables .....96
- G.2 Displaying Chromatograms from Mass Spectra .....98

## **Appendix H Editing and Adding Batch Files During Sequential Analysis**

---

<b>H.1</b>	Editing Batch Files .....	99
<b>H.2</b>	Adding Batch Files (Batch Queue) .....	101
H.2.1	Creating Batch Files to Add .....	101
H.2.2	Adding Batch Files .....	102

## **Appendix I Reducing the Carrier Gas Flow Rate After Sequential Analysis and Using the Ecology Mode**

---

<b>I.1</b>	Reducing the Carrier Gas Flow Rate After Sequential Analysis .....	105
I.1.1	Creating a Method File That Reduces the Carrier Gas Flow Rate .....	105
I.1.2	Creating Batch Files .....	106
<b>I.2</b>	Ecology Mode (This feature applies to QP2010 Ultra and QP2010 SE models.) .....	109
I.2.1	Setting the Mode Manually .....	109
I.2.2	Setting the Mode Using Batch Processing .....	110

## **Appendix J Printing Reports**

---

<b>J.1</b>	Printing Images (Printing Spectra and Chromatograms Displayed in Windows) .....	112
<b>J.2</b>	Creating Reports .....	114
J.2.1	Using Templates .....	114
J.2.2	Using Previously Created Report Files .....	115
J.2.3	Manually Setting Report Content .....	116

# 1





# Overall Configuration of GCMSSolution

## 1.1 Programs

## 1


GCMSSolution is made up of the programs described below.

Select the program that is appropriate for the purpose (e.g., analysis or data processing).

Icon	Name	Description
 GCMS Real Time Analysis	GCMS Real Time Analysis	Used to start up and shut down the instrument, make configuration settings, and perform analysis.
 GCMS Analysis Editor	GCMS Analysis Editor	Used to create and edit method files and batch files during analysis.
 GCMS Postrun Analysis	GCMS Postrun Analysis	Used to perform qualitative and quantitative processing, print reports, and perform other tasks involving data processing.
 GCMS Browser	GCMS Browser	Used to perform qualitative and quantitative processing, print reports, and perform other data processing tasks for multiple data files.

 Ver. 2.5


This manual describes operating instructions primarily for new features added in GCMSSolution Version 2.6.

However, items indicated with " Ver. 2.5" can also be performed using previous methods. Use whichever procedure is easier.

 **NOTE**

Icons and windows for functions that can only be used on QP2010 Ultra or QP2010 SE model will not be displayed on the software if the GCMS model used is QP2010, QP2010 Plus or QP2010S.

## 1.2 Flowchart of Operating Procedure





# 2


## Starting GC/MS

### 2.1 Turning ON the Power


Switch ON any peripheral or accessory equipment connected to the system, before switching ON the main GC/MS system.

2

#### 1 Turn ON the power to the GC.



#### 2 Turn ON the power to the MS.




#### 3 Turn ON the power to the PC, printer, and display.




**4** Double-click the **GCMS Real Time Analysis** (GCMS Real Time Analysis) icon.  
The [GCMS Real Time Analysis] program starts.


**5** Click [OK].



## 2.2 Layout of Operating Areas



**GCMS Real Time Analysis**



**GCMS Postrun Analysis**

No.	Name	GCMS Program	Explanation
1	Title Bar	Real Time Analysis, Analysis Editor, Postrun Analysis, and Browser	Displays the name of the program, process, and file currently running or being processed.
2	Menu Bar	Real Time Analysis, Analysis Editor, Postrun Analysis, and Browser	Displays command menus corresponding to the window currently open.
3	Toolbar	Real Time Analysis, Analysis Editor, Postrun Analysis, and Browser	Displays command tool buttons corresponding to the window currently open.

No.	Name	GCMS Program	Explanation
4	Assistant Bar	Real Time Analysis, Analysis Editor, Postrun Analysis, and Browser	Command icons are arranged in order of typical operation sequence. The assistant bar is named according to the window that is currently open. For example, when the [Batch] window is open, the assistant bar is named the [Batch] assistant bar.
5	Instrument Monitor	Real Time Analysis	Displays analytical instrument parameter values in real time.
6	Data Explorer	Real Time Analysis, Analysis Editor, Postrun Analysis, and Browser	Used to easily load analytical data or method files. It lists files in the selected folder, according to file type.

**NOTE**

The assistant bar, instrument monitor, and Data Explorer can be shown or hidden by selecting [Show/Hide] on the [View] menu.


2

## 2.3 Inspecting Consumable Items and Maintenance Parts

Check the state of the GC/MS consumable items using the procedure described below.

- 1** Move the mouse pointer over the icon for a consumable item in the instrument monitor to display the current state and the recommended replacement point for the corresponding item.

When a consumable item approaches its recommended maximum usage frequency, the background of the corresponding icon turns black to alert the user.



This note is shown when mouse pointer is moved over the septum icon. This means that the septum has been used 13 times out of a maximum 100 times.

**NOTE**

When replacing the analysis column, or when a consumable item has passed its recommended replacement point, perform maintenance with reference to ["Appendix D Maintenance" P.83](#).

Depending on the analysis content, the appropriate replacement frequency may be greater than the recommended frequency.

## 2.4 System Configuration


Check and set the modules used for analysis using the procedures described below.

### 2.4.1 Setting the Modules Used for Analysis


- 1 Click the [System Configuration] icon on the [Real Time] assistant bar.  
The [System Configuration] window opens.



- 2 Check that the components shown in the [Modules Used for Analysis] area correspond to the actual modules in the GC/MS system that are to be used for the analysis.





If the modules to be used for current analysis do not correspond to the modules shown in this window, set as shown in the following example:

- 1 Select [AOC-20i+s] in the [Available Modules] area if for example, AOC-20i with AOC-20s are to be used for analysis.
- 2 Click  to register the module in [Modules Used for Analysis].

## 2.4.2 Checking Column Information


Information for the column attached to the instrument must be specified correctly. Check or change the column information settings after installing a column.

- 1 Double-click the  (Column) icon in [Modules Used for Analysis]. The [Modules of Analytical Line#] window opens.




If the column to be used is not displayed under [Selected Column] area:

- 1 Click the name of the column to be used in the [Registered Columns] table.



- 2 Click [Select].  
The column information is displayed under [Selected Column].



**If the column to be used is not registered:**

- 1 Click [Add].  
A row is added.

Modules of Analytical Line#1

AOC-20i | SPL1 | Column | MS

Please input the column information on the [Registered Columns] table, and click the [Select] button.

Selected Column

Name : Rtx-5MS  
 Serial # : 904661      Thickness : 0.25 um  
 Length : 30 m      Diameter : 0.25 mm  
 Max. Usable Temp. : 350 °C      Installation Date : 5/11/2009

Registered Columns

	Column Name	Serial #	Thickness (um)	Length (m)	Diameter (mm)
1	Rtx-624	657772	1.8	60	0.
2	Rtx-5MS	904661	0.25	30	0.
3			0	0	

Description :

Select Add Delete

- 2 Enter the column information.  
Column information is usually indicated on the packaging box of the column, on the column tag, or in the column specification sheet that is usually included inside the column box.

Modules of Analytical Line#1

AOC-20i | SPL1 | Column | MS

Please input the column information on the [Registered Columns] table, and click the [Select] button.

Selected Column

Name : Rtx-5MS  
 Serial # : 904661      Thickness : 0.25 um  
 Length : 30 m      Diameter : 0.25 mm  
 Max. Usable Temp. : 350 °C      Installation Date : 5/11/2009

Registered Columns

	Column Name	Serial #	Thickness (um)	Length (m)	Diameter (mm)
1	Rtx-624	657772	1.8	60	0.
2	Rtx-5MS	904661	0.25	30	0.
3	Rxi-5SII MS		0	0	


Description :

Select Add Delete

**NOTE**


The maximum temperature that can be specified for the column oven, interface, and sample injection unit is normally the [Max. Usable Temp] setting entered here. For the protection of columns, always enter a [Max. Usable Temp] setting.

- 2** Click [OK].  
The [System Configuration] window returns.



### 2.4.3 Enabling the Modules Used for Analysis

- 1** Click [Set].  
The system configuration information is transferred to the instrument.




## 2.5 Vacuum System Startup


Open the carrier gas cylinder valve to supply carrier gas.

If carrier gas is being controlled by accessory/peripheral equipment, use that equipment to supply carrier gas before starting the vacuum system.


- 1 Click the [Vacuum Control] icon on the [Real Time] assistant bar.  
The [Vacuum Control] window opens.



- 2 Click [Auto Startup].  
The vacuum system starts.




- 3 When [Completed] is displayed, click [Close].





## 2.6 Checking for Vacuum Leakage

- 1** Wait for 10 minutes after starting up the vacuum system.
- 2** Click the [Tuning] icon on the [Real Time] assistant bar.  
The [Tuning] window opens.




2

- 3** Click the [Peak Monitor View] icon on the [Tuning] assistant bar.  
The [Peak Monitor] window opens.



## 4 Check for vacuum leaks.

The screenshot shows the GCMS Real Time Analysis software interface. The main window displays three mass spectra plots for m/z 18, 28, and 32. The m/z 28 peak is significantly higher than the m/z 18 peak. The interface includes a 'Monitor Group' dropdown menu, a 'Detector' voltage control, and various system parameters like 'Low Vacuum' and 'High Vacuum'. The 'Detector' voltage is set to 0.65 kV. The 'Low Vacuum' is 5.6e+000 Pa and the 'High Vacuum' is 1.4e-004 Pa. The 'Detector' voltage is set to 0.65 kV. The 'Low Vacuum' is 5.6e+000 Pa and the 'High Vacuum' is 1.4e-004 Pa.

- 1 Click the arrow button in [Monitor Group] setting, and select [Water, Air] from the list.
- 2 Click  (Filament ON/OFF) to turn ON the filament. Peaks will be displayed in the three windows.
- 3 Change the detector voltage gradually by clicking the up or down arrow buttons so that the peak height for  $m/z$  18 (water) corresponds to half the height of the display window.
- 4 Compare the peak heights for  $m/z$  18 (water) and  $m/z$  28 (nitrogen).  
Check that the peak height for  $m/z$  28 (nitrogen) is not more than twice that for  $m/z$  18 (water).

### NOTE

If the peak height for  $m/z$  28 (nitrogen) is more than twice that for  $m/z$  18 (water), it is possible that there is an air leak. Search for the location of the leak.


### Reference

Refer to the System User's Guide for details on how to check for vacuum leaks.

- 5 Click  (Filament ON/OFF) to turn OFF the filament.

## 5 Close the [Tuning] window.

The message [Save current tuning file?] is displayed. Click [No].



## 2.7 Autotuning

Wait for approximately 2 hours (before starting qualitative analysis) or 4 hours (before starting quantitative analysis) after starting up the vacuum system and then perform autotuning using the procedures described below.

### 2.7.1 Setting Analysis Conditions

If no analysis conditions have been created, start from ["2.7.2 Executing Autotuning" P.14.](#)

If a method file is already created, parameters can be specified in the instrument according to the following procedure.


#### NOTE

However, parameters for an accessory or peripheral equipment, except for AOC-20 auto-injector/auto-sampler, cannot be specified by using the following procedure. When using an accessory/peripheral equipment, set the parameters on the equipment itself, or by using the software specific to that equipment/device.


- 1 Click the [Data Acquisition] icon on the [Real Time] assistant bar.  
The [Acquisition] window opens.




- 2 Click  (Open) on the toolbar.



- 3** Select the method file to load, then click [Open].  
The method file is loaded.



- 4** Select [Download Initial Parameters] on the [Acquisition] menu.




**NOTE**

If the message "The hardware configuration for this method is different from the current instrument configuration. The measurement condition in the method file is modified according to the current instrument configuration." appears, click [OK].

## 2.7.2 Executing Autotuning

- 1** Click the [Tuning] icon on the [Real Time] assistant bar.  
The [Tuning] window opens.




- 2** Click the **[Peak Monitor View]** icon on the **[Tuning]** assistant bar.  
The **[Peak Monitor]** window opens.



2


- 3** Select **[New Tuning File]** on the **[File]** menu.




- 4** Select **Tuning Mode** appropriate for the application.  
(This feature applies to QP2010 Ultra and QP2010 SE models.)

When creating a new tuning file, choose the tuning mode appropriate for the concentration level of target compounds being measured. Since the tuning file is created with an emission current corresponding to the selected mode, it enables measuring samples with an appropriate dynamic range.

- QP2010 Ultra: High concentration (20  $\mu$ A), standard (60  $\mu$ A, default), or high sensitivity (150  $\mu$ A)
- QP2010 SE: High concentration (20  $\mu$ A) or standard (60  $\mu$ A, default)



- 5** Select the filament to be used.




## 6 Click the [Start Auto Tuning] icon on the [Tuning] assistant bar.



## 7 Enter a file name and click [Save] to start autotuning.

When autotuning is completed, a report is printed.



## 8 Close the [Tuning] window.


The message [Save current tuning file?] is displayed. Click [Yes].

## 2.7.3 Checking Autotuning Results

### 1 Check the results of autotuning.

C:\GCMSolution\System\Tune1\20091115.qg

### Tuning Condition ###  
 Resolution Adjustment : YES  
 FWHM of Peak Profile : 0.6  
 Sensitivity Adjustment : YES  
 Target Mass Adjustment : 264.00  
 Mass Calibration : YES  
 Mass Pattern Adjustment : NO




MS : GCMS-QP2010 Ultra  
 Ionization Mode : EI  
 Tuning Date : 11/18/2009 10:49:35 AM  
 Filament# : 1  
 Lens1 : -1.8 V  
 Lens2 : -17.0 V  
 Lens3 : -2.1 V  
 Lens4 : -50.0 V  
 RF Gain : 5061  
 RF Offset : 4952

Detector : 0.86 kV

Ionization voltage : 70 V  
 Emission current : 60  $\mu$ A  
 Main-rod : -3.5 V  
 Conversion dynode : -10 kV

IonSourceTemp : 200 °C  
 Low Vacuum : 5.5e+000 Pa  
 High Vacuum : 1.3e-004 Pa  
 Interface Temp. : 250 °C  
 Oven Temp. : 45 °C  
 Column Pressure : 68 kPa  
 Column Flow : 1.2 mL/min  
 Column Diameter : 0.25 mm  
 Column Length : 30.0 m



- 1 Check that the FWHM (full width at half maximum) values are in the range 0.5 to 0.7.
- 2 Check that the detector voltage does not exceed 2 kV.
- 3 Check that the base peak values are 18 or 69.
- 4 Check that the relative intensity ratio for  $m/z$  502 is at least 2 % (for QP2010S and SE : 1 %).
- 5 Check that the peak intensity for  $m/z$  69 is at least twice that for  $m/z$  28.

#### NOTE

If any irregularities are discovered above, possible causes could include a vacuum leak, poor column connections, or contaminated ion source.

See "[Appendix D Maintenance](#)" P.83 to implement corrective measures.

# 3


## Creating a Folder


### 3.1 Creating a Folder with Data Explorer


- 1 Click the [Data Acquisition] icon on the [Real Time] assistant bar. The [Acquisition] window opens.



- 2 Click  (Data Explorer) on the toolbar to display Data Explorer.




- 3 Click  (Project (Folder) Selection). The [Project (Folder) Selection] window opens.






- 4** Click the **GCMSsolution** folder in the **C** drive in **Data Explorer** and click **[New Folder]**.  
The **[Create New Folder]** window opens.



3

- 5** Type a folder name and click **[OK]**.  
A folder is created in the **GCMSsolution** folder in the **C** drive and the **[Project (Folder) Selection]** window returns.



- 6** Click **[Close]**.  
To select an existing folder, see ["Appendix C Using Data Explorer" P.81](#).

# 4

# Qualitative Analysis

## 4.1 Creating a Method File

Set the instrument (i.e., autosampler, GC, MS) parameters and similarity search parameters using the procedure described below.



### NOTE

Use default values for parameters that are not covered by the following explanations.

1


Click the **[Data Acquisition]** icon on the **[Real Time]** assistant bar.

The **[Acquisition]** window opens.




2

Select **[New Method File]** on the **[File]** menu.



## 4.1.1 Setting Autosampler Parameters

- 1 Click the [Sampler] tab and specify the number of rinses appropriate for the sample.



4

## 4.1.2 Setting GC Parameters

- 1 Click the [GC] tab and set the analysis conditions.

	Rate	Final Temperature	Hold Time
0	-	45.0	1.00
1	45.00	130.0	0.00
2	12.00	180.0	0.00
3	7.00	240.0	0.00

- 1 Input an initial temperature for the column oven (40 to 100 °C).
- 2 Input an injection temperature based on consideration of the boiling point of the target compound (200 to 300 °C).
- 3 Select [Split] or [Splitless].

**NOTE**

## Selecting Injection Mode

- Split: Select this mode if the concentration of the target compound is high. As a rough guideline, select this mode when the target compound concentration is greater than 10 ng/uL.
- Splitless: Select this mode if the concentration of the target compound is low. As a rough guideline, select this mode when the target compound concentration is less than 10 ng/uL.

4 Select [Pressure] when the method calls for a constant pressure mode, and select [Linear Velocity] when the method calls for a constant linear velocity mode for the carrier gas. When no reference method is available, select [Linear Velocity].

5 When no reference method is available, refer to the table "Typical Pressure Settings for Carrier Gas" to set an initial value for the pressure. The linear velocity will be set automatically.

## Typical Pressure Settings for Carrier Gas

Middle bore capillary column (I.D. 0.25 mm)		Semi-wide bore capillary column (I.D. 0.32 mm)	
30 m	60 m	30 m	60 m
75 to 150 kPa	100 to 250 kPa	30 to 50 kPa	50 to 100 kPa

6 If "Split" is selected as the injection mode, enter a split ratio. If "Splitless" is selected, enter "-1.0".

7 Set appropriate conditions for separating the target compound from other peaks.

## 4.1.3 Setting MS Parameters

# 1

Click the [MS] tab and set the analysis conditions.

The screenshot shows the MS parameter configuration window. Key elements include:

- MS Tab:** Selected at the top.
- Interface Temp.:** 250 °C (Callout 1).
- Solvent Cut Time:** 2.5 min (Callout 2).
- GC Program Time:** 26.29 min (Callout 3).
- Acq. Mode:** Scan (Callout 5).
- Start m/z:** 35.00 (Callout 6).
- End m/z:** 500.00 (Callout 6).
- Table:**

Group#1 - Event#1	Start Time (min)	End Time (min)	Acq. Mode	Event Time(sec)	Scan Speed	Start m/z	End m/z	Ch1 m/z	Ch2 m/z	Ch3 m/z	Ch4 m/z
1	3.00	26.00	Scan	0.30	1666	35.00	500.00				
2	0.00	0.00	Scan	0.00	0	0.00	0.00				

1 Input [Interface Temp.] (200 to 300 °C).


2 Input [Start Time] and [End Time] according to the note below.

**NOTE**

In the absence of information about the elution time of the solvent peak, set [Start Time] to zero minutes, and set [End Time] to the [GC Program Time] value. After one analysis of a standard sample or the solvent, and obtaining the solvent peak profile, change the [Start Time] to a time after the end of the solvent peak (see the figure shown on page 23).

- 3 Click [Relative to the Tuning Result].  
If peak intensity is too low, change the value within the range +0.1 to +0.3., as necessary.
- 4 Input a value that is 0.5 minutes less than the [Start Time] setting. (If the resulting value is less than zero, enter "0".)

Relationship between Start Time and Solvent Elution Time




- 5 Select [Scan].
- 6 Enter the mass range to be measured, where [Start m/z] is the lower mass limit, and [End m/z] is the upper mass limit. The typical value for [Start m/z] is 35, and the typical value of [End m/z] is the highest molecular weight of the target compounds in the sample plus some margin of error (+15).


## 4.1.4 Setting Similarity Search Parameters

### 1 Select [Qualitative Parameters] on the [Method] menu.


The [Qualitative Parameters] window opens.



## 2 Click the [Similarity Search] tab and set the search conditions.




- 1 Click  .  
The [Open File] window opens.




- 2 Open the library to be used.
- 2 To remove a library from the selection, highlight the library file name by dragging the mouse over it, then press the [Delete] key.
- 3 Select [Do not include duplicate hits].
- 4 After completing the settings, click [OK] to return to the original window.

## 4.1.5 Saving the Method File

- 1 Select [Save Method File As] on the [File] menu.



- 2 Enter a file name and click [Save].



---


## 4.2 Repeating Autotuning

If autotuning has not been performed under the analysis conditions, perform the procedures described under "[2.7 Autotuning](#)" P.13.

## 4.3 Sequential Analysis

Create a batch file necessary for qualitative analysis and perform sequential analysis using the procedures described below.


### 4.3.1 Creating a Batch File

To use an existing batch file, follow the procedure starting in "[4.3.2 Editing a Batch File](#)  Ver. 2.5" P.29.

- 1 Click the [Batch Processing] icon on the [Real Time] assistant bar.  
The [Batch Table] window opens.



- 2 Select [New Batch File] on the [File] menu.







### 3 Click the [Wizard] icon on the [Batch] assistant bar.


The [Batch Table Wizard] window opens.




### 4 With the Batch Table Wizard, make the appropriate settings and create a batch table.



- 1 Click [New].
- 2 Click [Unknown Only].
- 3 Click  and specify the method file to be used.
- 4 Deselect both [Data Processing] items.
- 5 Click [Next].



- 6 Input [Vial #] and [Sample Count].
- 7 Input [Injection Volume].
- 8 Click [Next].




- 9 Enter [Data File Name].  
If the file name ends with a number, the files are named sequentially.
- 10 Click [Finish]. The batch table is displayed.

## 4.3.2 Editing a Batch File Ver. 2.5

For routine analyses, it may be more convenient to partially edit and execute existing batch files. The following procedures describe how to edit information in specified row(s) collectively.

### 1 Double-click the batch file to be edited.



### 2 Add or delete rows depending on the number of samples being analyzed.

Folder: C:\GCMSsolution\20091115_operation guide							
	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	1	Methylene Chloride		0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	2	PAH 0.005ppm	STD-0001	1:Standard(I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	3	PAH 0.01ppm	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	4	PAH 0.05ppm	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	5	PAH 0.1ppm	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	1	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample01.qgd
7	6	river water 1	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample02.qgd
8	7	river water 2	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample03.qgd

- 1 Click on the row number to be edited to highlight the whole row. Multiple rows can be selected by dragging the mouse over multiple row numbers.

Folder: C:\GCMSsolution\20091115_operation guide							
	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	1	Methylene Chloride		0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	2	PAH 0.005ppm	STD-0001	1:Standard(I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	3	PAH 0.01ppm	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	4	PAH 0.05ppm	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	5	PAH 0.1ppm	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	1	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample01.qgd
7	6	river water 1	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample02.qgd
8	7	river water 2	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample03.qgd

- 2 Right-click on the selected row, and select the appropriate editing command from the menu that is displayed.

Copy Row  
Add Row  
Insert Row  
Paste Row  
Delete Row

Menus	Explanation
Copy Row	Copies the selected row.
Add Row	Adds a row to the end.
Insert Row	Inserts a new row above the selected row.
Paste Row	Pastes the copied row.
Delete Row	Deletes the selected row.


In addition, rows can be added easily by repeatedly entering "1" as the vial number and pressing the down arrow key on the keyboard.

The following steps can be performed to quickly copy the contents of one row onto several rows, or to fill in the contents of several rows quickly. Serial numbers will be added to the end of the sample name and data file name, for examples, in the edited rows.

### 3 Edit the vial number, sample name and data file name in an existing row, which is to be the first of several rows to be edited.

Folder: C:\GCMSsolution\20091115_operation guide							
	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	1	Methylene Chloride		0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	2	PAH 0.005ppm	STD-0001	1:Standard(I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	3	PAH 0.01ppm	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	4	PAH 0.05ppm	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	5	PAH 0.1ppm	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	1	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample01.qgd
7	6	river water 1	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample02.qgd
8	7	river water 2	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample03.qgd
9	1			0:Unknown	IT QT		


- 1 Directly enter the vial number and sample name. Add a number at the end of sample name.
- 2 To set the data file name, click within the cell, then click the arrow buttons that appear. The [Select Data File] window opens. Directly enter the data file name in the [Select Data File] window, then click [Open].



#### NOTE

It is not necessary to type file extensions when entering data file names.

- 3 To set the method file, click within the cell, then click the arrow buttons that appear. The [Select Method File] window opens. Select the method file to use.



## 4 Drag the mouse from the edited row to the row specified with serial numbers.

Folder: C:\GCMSsolution\20091115\_operation guide

	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	1	Methylene Chloride		0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	2	PAH 0.005ppm	STD-0001	1:Standard(I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	3	PAH 0.01ppm	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	4	PAH 0.05ppm	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	5	PAH 0.1ppm	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	1	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample01.qgd
7	6	river water 1	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample02.qgd
8	7	river water 2	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample03.qgd
9	1			0:Unknown	IT QT		
10	1			0:Unknown	IT QT		
11	1			0:Unknown	IT QT		
12	1			0:Unknown	IT QT		
13	1			0:Unknown	IT QT		



### NOTE

To collectively edit specified cells, without changing other settings, drag the mouse across the cells to edit and perform the operation described in step 6.


For example, this can be used to edit the vial numbers and data file names for standard samples used to create calibration curves (for quantitative analysis).

Folder: C:\GCMSsolution\20091115\_operation guide

	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	Methylene Chloride		0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	PAH 0.005ppm	STD-0001	1:Standard(I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	PAH 0.01ppm	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	PAH 0.05ppm	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	PAH 0.1ppm	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample01.qgd


## 5 Select [Fill Down] on the [Edit] menu.

The entire content of the first row is copied.




## 6 Select [Fill Series] on the [Edit] menu.

Edited parameters will be appended with serial numbers.




### 4.3.3 Saving Batch Files

- 1 Select [Save Batch File As] on the [File] menu.



- 2 Open the folder where the method file is saved, enter a name, and save the file.





## 4.3.4 Executing Sequential Analysis

- 1 Set the syringe rinse solvent and samples in the autosampler.
- 2 Click the [Start] icon on the [Batch] assistant bar.  
Analysis starts.




### NOTE

- To abort batch processing, click the  (Stop) icon on the [Batch] assistant bar.
- To modify or add batch files while analysis is in progress, see ["Appendix H Editing and Adding Batch Files During Sequential Analysis" P.99](#).
- To execute only specified rows, select the rows by clicking or dragging the mouse, then start the analysis.



## 4.4 Analyzing Data

Use the procedure described below to perform basic qualitative data processing for data measured in Scan mode, for examples, to display mass spectra, perform background subtraction, and perform similarity search.


- 1** Double-click the  (GCMS Postrun Analysis) icon.  
The [GCMS Postrun Analysis] program starts.

- 2** Click the [Qualitative] icon on the [Postrun] assistant bar.



### 4.4.1 Loading Data Files

- 1** With reference to "[Appendix C Using Data Explorer](#)" P.81, double-click the data file to be analyzed.  
The data file to be analyzed opens.






## 4.4.2 Zoom in on a Peak in the Chromatogram and Displaying Mass Spectra

- 1 Specify a range in the TIC window by dragging the mouse so that both the peak top and baseline are highlighted.

Drag the mouse so that both the peak top and baseline are displayed.





The top window displays the full chromatogram, and the second window displays any enlarged portion of the same chromatogram.


### NOTE

To undo the zoom, right-click in the MC window and select [Undo Zoom] on the pop-up menu.


- 2** Move the mouse pointer to the peak top and double-click.  
The mass spectrum for the peak top is displayed.



If a mass spectrum location other than the peak top is displayed, align the bar position with the peak top by clicking .




#### NOTE


If a red peak appears in the mass spectrum, it indicates that the signal intensity for this m/z has exceeded the scale. To display the proper mass spectrum for target compounds, click the left or right arrow buttons  to select a mass spectrum that shows no red peaks.

### 4.4.3 Removing the Background

The mass spectrum corresponding to a peak in the chromatogram has a contribution from the background. This background mass spectrum may interfere with the subsequent qualitative analysis, hence it needs to be subtracted from the compound's mass spectrum. The steps shown in this section describe the procedure for removing the background mass spectrum from a compound's mass spectrum.


- 1 Click  (Spectrum Subtraction) on the toolbar and position the mouse pointer in the MC window.

A bar is displayed.



- 2 Double-click at the background processing position (see NOTE below).


The target spectrum, from which the background spectrum has been subtracted, is displayed.




#### NOTE

With the following types of peaks, process the parts indicated by arrows as background.


Example 1)



Example 2)



Example 3)





Background spectrum can be subtracted from one of positions.

## 4.4.4 Searching for Similarity Ver. 2.5

- 1 Click the [Similarity Search] icon on the [Qualitative] assistant bar. The [Similarity Search Results] window opens.



- 2 Click the  buttons on the left of the mass spectra as appropriate and check the mass spectra for the compounds found.



No.	Explanation
1	Similarity: The closer this value is to 100, the greater the similarity in mass spectra.
2	To enter a compound name in the spectrum table, select the box for the applicable compound.
3	Hit numbers for the compounds found.
4	Use to switch between the mass spectra for the compounds found.

### 3 After checking the mass spectra, click (Register Target Spectrum to Spectrum Process Table).

The mass spectrum is registered.



#### NOTE

Registering identification results in the spectrum process table allows referencing those results at a later time or outputting them as a report.

### 4 Close the [Similarity Search Results] window.


## 4.4.5 Displaying Mass Chromatograms

In the following situations, displaying mass chromatograms (MC) makes it easier to analyze the data. For instructions of how to display mass chromatograms, see ["Appendix G Displaying Chromatograms" P.96](#).

4

### ■ Confirming the Purity of Peaks


Displaying mass chromatograms can be used to check the presence of two or more overlapping peaks, or in other words, to check the purity of a peak in the chromatogram.



## ■ Looking for Target Compound Peaks Among Multiple Peaks

In some cases, peaks for target compounds cannot be confirmed in a total ion chromatogram (TIC).

If characteristic mass spectral peaks (i.e.,  $m/z$ ) of the target compounds are known, displaying the mass chromatograms makes it easier to check the position of the target compound's peaks in the chromatogram.



### NOTE

If characteristic spectral peaks are not known for the target compounds, see ["Appendix F Index Searches" P.94](#) to check the mass spectra of target compounds.

## 4.4.6 Registering Spectra Displayed for Target Compounds

Ver. 2.5

- 1 Register all the other target compounds using the procedure described in ["4.4.2 Zoom in on a Peak in the Chromatogram and Displaying Mass Spectra" P.35](#).

## 4.4.7 Editing the Spectrum Process Table


- 1 Click the [Qualitative Table] icon on the [Qualitative] assistant bar. The [Qualitative Table] window opens.




2 Click  (Maximize).

Qualitative Table									
Spectrum					Background				
	Ret. Time	Start Tm	End Tm	Ret. Time	StartRT	EndRT	Search	Report	Name
1	...	4.860	4.870	4.865	4.835	4.315	Done	<input checked="" type="checkbox"/>	Naphthalene-D8
2	...	6.985	6.995	6.990	6.950	7.045	Done	<input checked="" type="checkbox"/>	Acenaphthylene
3	...	7.235	7.245	7.240	7.200	7.285	Done	<input checked="" type="checkbox"/>	Acenaphthene-d10
4	...	8.200	8.210	8.205	8.160	8.250	Done	<input checked="" type="checkbox"/>	Fluorene
5	...	10.250	10.260	10.255	10.205	10.280	Done	<input checked="" type="checkbox"/>	Phenanthrene-D10
6	...	10.305	10.315	10.310	10.280	10.365	Done	<input checked="" type="checkbox"/>	Anthracene
7	...	10.410	10.420	10.415	10.365	10.470	Done	<input checked="" type="checkbox"/>	Anthracene
8	...	14.210	14.220	14.215	14.160	14.290	Done	<input checked="" type="checkbox"/>	Fluoranthene
9	...	17.785	17.795	17.790	17.740	17.810	Done	<input checked="" type="checkbox"/>	Triphenylene
10	...	17.795	17.835	17.830	17.810	17.865	Done	<input checked="" type="checkbox"/>	Chrysene-D12
11	...	17.890	17.900	17.895	17.865	17.965	Done	<input checked="" type="checkbox"/>	Benzo[a]anthracene
12	...	20.290	20.300	20.295	20.250	20.320	Done	<input checked="" type="checkbox"/>	Benzo[e]pyrene
13	...	20.340	20.350	20.345	20.320	20.455	Done	<input checked="" type="checkbox"/>	Benzo[e]acephenanthylene
14	...	20.915	20.925	20.920	20.875	20.995	Done	<input checked="" type="checkbox"/>	Benzo[j]acephenanthylene
15	...	21.035	21.045	21.040	20.995	21.140	Done	<input checked="" type="checkbox"/>	Pyrene-D12
16	...	22.875	22.885	22.880	22.830	22.895	Done	<input checked="" type="checkbox"/>	Benzo[ghi]perylene
17	...	22.915	22.925	22.920	22.895	23.020	Done	<input checked="" type="checkbox"/>	p-Bis[phenylethynyl]benzene
18	...	23.325	23.335	23.330	23.290	23.410	Done	<input checked="" type="checkbox"/>	Benzo[ghi]perylene

3 Double-click the first row in the spectrum table.  
The [Similarity Search Results] window opens.

4 Click the  buttons on the left of the mass spectra as appropriate to double-check the mass spectra for the compounds found.




No.	Explanation
1	Similarity: The closer this value is to 100, the greater the similarity in mass spectra.
2	To enter a compound name in the spectrum table, select the box for the applicable compound.
3	Click to copy the selected compound names to the spectrum table.

No.	Explanation
④	Use to switch between the mass spectra for the compounds found.
⑤	Hit numbers for the compounds found.
⑥	Allows switching between search results for each row in the spectrum process table.


**5** After checking the mass spectra, close the [Similarity Search Results] window.


**6** To edit compound names or other information in the spectrum process table, click the desired row, then select the [Edit] menu.



**7** Close the [Qualitative Table] window.

#### 4.4.8 Saving Data Files

**1** Click  (Save) on the toolbar.  
The qualitative table is saved in the data file.






## 4.5 Printing Qualitative Analysis Reports


The procedures are described below on how to use a template to create a report of analyzed data, how to edit the area of the chromatogram to display in the report, and how to edit the number of compounds to be displayed in the report of similarity search results.

### 4.5.1 Loading Report Formats


- 1 Click the **[Report]** icon on the **[Qualitative]** assistant bar.  
The **[Data Report]** window opens.



- 2 Click **[New Format File]** on the **[File]** menu.  
The **[File New]** window opens.




- 3 Select **[Use Template]** and select the format **[Qualitative Analysis Report]**.




- 4 Click **[OK]**.  
The **[Qualitative Analysis Report]** format opens.

## 4.5.2 Editing Report Formats


- 1** Double-click on the chromatogram.  
The [GCMS Chromatogram Properties] window opens.



- 2** Click the [Chromato] tab.




### 3 In the [Area] area, deselect [Auto] for the X-axis and enter the time range.




4

### 4 Click [OK].


### 5 Click the next page icon on the toolbar to display the second page.




### 6 Double-click on the [Library] display item. The [GCMS Library Properties] window opens.



## 7 Click the [Result] tab.




## 8 Enter the [Maximum Compound Number] (maximum number of search results to display).



## 9 Click [OK].


## 4.5.3 Outputting Reports


- 1 Click the [Preview] icon on the [Data Report] assistant bar.  
The print preview window opens.



- 2 After checking the report content, click [Print] to print the report.

The screenshot shows the GCMS Postrun Analysis (Admin) window. The title bar reads 'GCMS Postrun Analysis (Admin) [Data Report - ...004.qgd(Report in Data File)]'. The window contains two side-by-side report pages. The left page is titled 'Qualitative Analysis Report' and dated '11/17/2009'. It displays 'Sample Information' and a chromatogram. The right page is also titled 'Qualitative Analysis Report' and dated '11/17/2009'. It displays a list of peaks with their retention times and chemical structures. The window has a standard toolbar with 'Print', 'Next Page', 'Prev Page', 'One Page', 'Zoom In', 'Zoom Out', and 'Close' buttons. The status bar at the bottom indicates 'Pages 1-2 / Total 19' and 'NUM'.

- 3 Click  (Save) on the toolbar.  
The report is saved as part of the data file.



# 5

# Quantitative Analysis

## 5.1 Creating a Method File

With reference to "4 *Qualitative Analysis*" P.20, analyze standard samples and register the retention times and mass spectra of the target compounds in the spectrum process table.


### 5.1.1 Creating a Compound Table

- 1 Start the [GCMS Postrun Analysis] program and click the [Compound Table] icon on the [Postrun] assistant bar.

The [Create Compound Table] window opens.




- 2 From Data Explorer, double-click the data file in which the spectrum process table for the target compounds was saved.




- 3 Click the [Wizard (New)] icon on the [Compound Table] assistant bar. The [Compound Table Wizard] window opens.




## 4 Select [Use current Spectrum Process Table] and then click [Next].




## 5 Click [Next].



## 6 Select a row in the table, check the mass spectrum for each compound, and click [Next].



## 7 Specify the calibration curve type, the quantitative method, and other parameters as required, and click [Next].



No.	Item	Explanation
1	Quantitative Method	<ul style="list-style-type: none"> <li>External Standard: Quantitation is performed using a calibration curve obtained from the absolute quantity (concentration) and the area or height value of the target compound in a standard sample.</li> <li>Internal Standard: An internal standard is added to the sample, the sample is analyzed, and quantitation is performed using the relationship between the relative sensitivity and the quantitative ratio with respect to the internal standard compound.</li> </ul>
2	Calculated by	Select [Area] or [Height]. Normally, select [Area].
3	# of Calib. Levels	Input the number of concentration levels of the calibration curve.
4	Unit	Set the concentration unit used for reports.
5	Format of Concentration	Set the number of digits used to indicate concentrations.



## 8 Make the appropriate settings for concentrations and measurement ions, and click [Next].

No.	Item	Explanation
1	Standard	Set the concentrations of the standard samples. If the concentration varies with the compound, make the necessary corrections after completing the wizard procedure.
2	Internal Standard	Set the concentration of the internal standard.
3	# of Reference Ions	Input the number of reference ions used to perform peak identification.
4	Decimal for mass	Determine the number of decimal places for target ion and reference ion m/z values. Selecting [1 Decimal] increases the sensitivity level.


5

## 9 Set the type, compound name, target ion, and reference ion for each substance. After entering the required information for all the compounds, click [Next].

- 1 Select [Target] in the [Type] list.  
Select [I.S.] when setting for an internal standard.
- 2 Change the type and m/z value.

- To change the type, click cell for the type to be changed and select "Target Ion", "Ref. Ion", or "Not used".
- To change the  $m/z$  value, click the cell for the  $m/z$  value to be changed and click the arrow button that appears.

The mass spectrum is displayed. To set the  $m/z$  value, enlarge the area near the spectrum peak to be set by dragging the mouse and double-click on the spectrum peak. Switching the compound ID number also switches the mass spectrum.



- 3 The names registered in the spectrum table are entered automatically. If there are blank spaces, enter the compound names.
- 4 Change the compound displayed by changing the ID number.

## 10 Click [Finish].

A compound table is created. Correct the contents of the compound table as required.

ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit	Ref.Ions	Conc.1	Conc.2	Conc.3	Conc.4	Event	STD Spec
1	Naphthalene	ISTD	1	136.00	4.865	0	mg/L	108.00-13	0.1	0.1	0.1	0.1	1	Registered
2	Acenaphthyle	Target	1	152.00	6.990	0	mg/L	151.00-15	0.005	0.01	0.05	0.1	1	Registered
3	Acenaphthen	ISTD	2	164.00	7.240	0	mg/L	162.00-16	0.1	0.1	0.1	0.1	1	Registered
4	Fluorene	Target	2	166.00	8.200	0	mg/L	165.00-16	0.005	0.01	0.05	0.1	1	Registered
5	Phenanthrene	ISTD	3	188.00	10.250	0	mg/L	189.00-18	0.1	0.1	0.1	0.1	1	Registered
6	Phenanthrene	Target	3	178.00	10.305	0	mg/L	176.00-15	0.005	0.01	0.05	0.1	1	Registered
7	Anthracene	Target	3	178.00	10.410	0	mg/L	176.00-17	0.005	0.01	0.05	0.1	1	Registered
8	Pyrene	Target	4	202.00	14.215	0	mg/L	200.00-20	0.005	0.01	0.05	0.1	1	Registered
9	Benzo[a]anthr	Target	4	228.00	17.790	0	mg/L	226.00-22	0.005	0.01	0.05	0.1	1	Registered
10	Chrysene-D12	ISTD	4	240.00	17.830	0	mg/L	236.00-24	0.1	0.1	0.1	0.1	1	Registered
11	Chrysene	Target	4	228.00	17.895	0	mg/L	226.00-22	0.005	0.01	0.05	0.1	1	Registered
12	Benzo[b]fluor	Target	4	252.00	20.295	0	mg/L	250.00-25	0.005	0.01	0.05	0.1	1	Registered
13	Benzo[k]fluor	Target	5	252.00	20.340	0	mg/L	250.00-25	0.005	0.01	0.05	0.1	1	Registered
14	Benzo[a]pyre	Target	5	252.00	20.920	0	mg/L	250.00-25	0.005	0.01	0.05	0.1	1	Registered
15	Perylene-D12	ISTD	5	264.00	21.035	0	mg/L	260.00-26	0.1	0.1	0.1	0.1	1	Registered
16	Indene[1,2,3-	Target	5	276.00	22.880	0	mg/L	274.00-27	0.005	0.01	0.05	0.1	1	Registered
17	Dibenz[a,h]an	Target	5	278.00	22.915	0	mg/L	279.00-13	0.005	0.01	0.05	0.1	1	Registered
18	Benzo[ghi]per	Target	5	276.00	23.330	0	mg/L	274.00-27	0.005	0.01	0.05	0.1	1	Registered
19		Target	5	TIC	0.000	0	mg/L		0.005	0.01	0.05	0.1	1	Registered

### NOTE


To correct the compound table, enter edit mode by clicking Edit at the top-right corner of the table. When editing is completed, return to display mode by clicking View.

## 11 Click the [Save Compound Table] icon on the [Compound Table] assistant bar.

The method file that was used to acquire the data will be selected automatically.



## 12 Click [Save].



### NOTE

If greater sensitivity is required, use the following procedure to create a quantitative analysis method for the SIM mode.

This completes the procedure for creating a quantitative method for Scan mode.


5

## 5.1.2 Creating a SIM Table

- 1 Click the [Create SIM Table [COAST]] icon on the [Compound Table] assistant bar. The [Select Method File] window opens.




- 2 Enter a file name and click [Save]. The [Create SIM Table [COAST]] window opens.



- 3 Click  (Maximize) in the [Creation of Automatic SIM (Scan/SIM) Table [COAST]] window.

## 4 A SIM table is created automatically. Check the chromatogram and SIM table and, if necessary, modify the table with reference to the following procedure.



### Modification Procedure

To ensure sufficient sensitivity, it is best to specify no more than 20  $m/z$  values per row (i.e. per group). If necessary, modify the SIM table.

No.	Item	Operation
①	Start Time	Click within the frame next to [Start Time], then point the mouse pointer to the location on the chromatogram corresponding to the time when data acquisition is to be started, then click the left mouse button. Normally start time is set 0.5 to 1 minute before elution of the first target compound.
②	End Times and $m/z$ Values for Each Group	Click the target row (i.e. group) and drag the mouse on the chromatogram to specify and enlarge the desired area. Next, click near the center of peaks labeled with compound names to specify the measurement end time for that group. This automatically sets the corresponding $m/z$ values for the group at the same time.

### NOTE

- To edit table rows (i.e. groups), right-click on the desired row and select the following on the menu that appears.
  - Add Row : Adds a row to the bottom of the table.
  - Insert Row : Inserts a new row above the selected row.
  - Delete Row : Deletes the selected row.
- To undo enlarging the chromatogram, right-click on the chromatogram and select [Undo Zoom] on the menu that appears.
- To split groups, use the following procedure. (Example: Splitting Group 3 into two groups)
  - Click the third row of the SIM table.
  - Right-click on the table and select [Insert Row].
  - Click the inserted row and drag the mouse on the chromatogram to specify and enlarge the desired area.
  - Click near the center of peaks labeled with compound names. Group 3 is divided into two groups.

- 5** When finished, click [OK].  
A method is created for SIM mode quantitative analysis.

## 5.2 Sequential Analysis

Create a batch file necessary for quantitative analysis and perform sequential analysis using the procedure described below.

### 5.2.1 Creating a Batch File


To use an existing batch file, follow the procedure starting in ["5.2.2 Editing a Batch File" Ver. 2.5" P.58](#).

- 1** Start the [GCMS Real Time Analysis] program and click the [Batch Processing] icon on the [Real Time] assistant bar.

The [Batch Table] window opens.




- 2** Select [New Batch File] on the [File] menu.





- 3** Click the [Wizard] icon on the [Batch] assistant bar.  
The [Batch Table Wizard] window opens.




## 4 Make the appropriate settings with the Batch Table Wizard and thereby create a batch table.



- 1 Select [New].
- 2 Select [Standard & Unknown].
- 3 Click  and specify the method file to be used.
- 4 Select [Quantitative].
- 5 Click [Next].




- 6 Input [Vial #].  
The number of calibration points is loaded automatically from the method.
- 7 Input [Injection Volume].
- 8 Input [Average Count] (i.e., the number of repetitions).
- 9 Click [Next].




- 10 Enter [Data File Name].  
If the file name ends with a number, the files are named sequentially.
- 11 Click [Next].

5



- 12 Input [Vial #].
- 13 Input [Sample Count].
- 14 Input [Injection Volume].
- 15 Click [Next].




- 16 Enter [Data File Name].  
If the file name ends with a number, the files are named sequentially.
- 17 Click [Finish].  
The batch table is displayed.

## 5.2.2 Editing a Batch File Ver. 2.5

- 1 To edit a batch file, see "[4.3.2 Editing a Batch File !\[\]\(48a7667d09d5a06397e047ee4537bb6f\_img.jpg\) Ver. 2.5](#)" P.29.


## 5.2.3 Saving Batch Files

- 1 Select [Save Batch File As] on the [File] menu.





## 2 Open the folder where the method file is saved, enter a name, and save the file.




### 5.2.4 Executing Sequential Analysis


- 1 Set the syringe rinse solvent and samples in the autosampler.
- 2 Click the [Start] icon on the [Batch] assistant bar.  
Analysis starts.

5



#### NOTE

- To abort batch processing, click the  (Stop) icon on the [Batch] assistant bar.
- To modify or add batch files while analysis is in progress, see ["Appendix H Editing and Adding Batch Files During Sequential Analysis" P.99](#).
- To execute only specified rows, select the rows by clicking or dragging the mouse and start the analysis.




## 5.3 Analyzing Data


### 5.3.1 Checking and Correcting Calibration Curves

- 1 Start the [GCMS Postrun Analysis] program and click the [Calibration Curve] icon on the [Postrun] assistant bar.


The [Calibration Curve] window opens.



- 2 Double-click the method file used in analysis from Data Explorer.



- 3** Select a compound in the compound table and click the calibration curve level.  
Check the calibration curve created and the chromatogram.



### Reference


If no peaks are identified or detected, perform identification or peak integration with reference to ["Manual Identification and Manual Peak Integration" P.63](#).




To change the method used to plot calibration curves, see ["Changing Parameters for Quantitative Processing" P.65](#).

### NOTE


Peaks that are detected in the chromatograms after automatic peak integration, will have peak detection marks (  $\uparrow\downarrow$  ).


The detected peaks are subjected to identification based on the retention times and ion ratios (  $\nabla$  peak identification mark).



Chromatogram	Countermeasure
<p>No peaks are detected.</p> 	Perform manual peak integration.
<p>Peaks are detected but different peaks are identified.</p> 	Perform manual identification.
<p>Peaks are detected and identified but peak integration is not performed properly.</p> 	Perform manual peak integration.

4

Only after correcting the calibration curves, click  (Save) on the toolbar to save the method file.




## Manual Identification and Manual Peak Integration

If no peaks are identified or detected, perform identification or peak integration using the procedure described below.

### Manual Identification


#### 1 Right-click in a chromatogram and select [Manual Identification] from the displayed menu.

A bar is displayed.



#### 2 Click the top of the peak to be identified.


The peak is identified.




## Manual Peak Integration

**1** Right-click in a chromatogram and select [Manual Peak Integrate...] from the displayed menu.


A bar is displayed.



**2** Drag the mouse from the start point to the end point of the peak. The [Base Line] window opens.



- 3 Select the baseline and click [OK].**  
The peak is integrated and identified.

**NOTE**

The same process can be accomplished by performing the following operations on the chromatogram.


Process	Operation	Explanation
Manual Identification	[Shift] + [Ctrl] + right-click	Identifies integrated peaks.
Manual Peak Integration	[Shift] + right-click-drag	Connects start point and end point as baseline.
Manual Peak Integration	[Ctrl] + right-click-drag	Connects points with horizontal baseline.

5


## ■ Changing Parameters for Quantitative Processing

Change quantitative processing parameters as necessary.


- 1 Click the [Quantitative Parameters] icon on the [Calibration] assistant bar.**  
The [Quantitative Parameters] window opens.



- 2 Click the [Quantitative] tab.**



### 3 Change the [Curve Fit Type], [Zero], and [Weighted Regression] settings, as necessary.



No.	Name	Explanation
①	Curve Fit Type	<p>Specifies how to plot the calibration curve.</p> <ul style="list-style-type: none"> <li>Linear: Determines the calibration curve as a straight line from the obtained values.</li> <li>Point to point: Points are connected by a broken line. No formula is displayed for point to point calibration curves.</li> <li>Quadratic: Curve is fit to each point using a quadratic equation. This requires at least three points on the calibration curve. For two points or less, the curve is calculated as linear.</li> <li>Mean RF: First, it determines straight lines passing through the origin and each point. Then it finds the simple average of the slopes for each line. Consequently, the resulting calibration curve always passes through the origin.</li> </ul>
②	Zero	Select either [Not Forced] or [Force Through]. Normally, select [Not Forced].
③	Weighted Regression	<p>A typical least squares method of plotting calibration curves could result in a quantitation error that is larger the lower the concentration at the calibration point. In general, when the calibration curve has a large dynamic range (maximum concentration is at least 50 times higher than the minimum quantitation limit), formulas are weighted to reduce the weight of higher concentration points of the calibration curve. Typically, formulas are optimized by checking the correlation coefficient and contribution ratio.</p> <ul style="list-style-type: none"> <li>[1/C<sup>2</sup>]: Formulas are weighted by the inverse of the concentration value squared.</li> <li>[1/C]: Formulas are weighted by the inverse of the concentration value.</li> <li>[1/A<sup>2</sup>]: Formulas are weighted by the inverse of the area value squared (or height value when a height is specified for the data used).</li> <li>[1/A]: Formulas are weighted by the inverse of the area value (or height value when a height is specified for the data used).</li> </ul>

### 4 After finishing making changes, click [OK].

Calibration curves are corrected according to the changed parameters.




## 5.3.2 Re-quantifying after Correcting a Calibration Curve

After correcting a calibration curve, re-quantify the data for samples with unknown concentrations.

- 1 Click the **[Batch Processing]** icon on the **[Postrun]** assistant bar.  
The **[Batch Table]** window opens.





- 2 Select **[New Batch File]** on the **[File]** menu.



- 3 Click the **[Select Data File]** icon on the **[Batch]** assistant bar.  
The **[Select Data File]** window opens.



- 4** Click the data file for sample with unknown concentrations, for which re-quantification is to be performed and click  (Add).  
The data file is selected.



- 5** Click [OK].

- 6** A batch table is displayed. Assign a name to the batch file and save it.

Folder: C:\GCMSsolution\20091115_operation guide								
	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Level#	ISTD Amt.
1	Unknown Sample	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_Sample1.qgd	1	1 1 1 1 1 1 1
2	Unknown Sample	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_Sample2.qgd	1	1 1 1 1 1 1 1
3	Unknown Sample	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_Sample3.qgd	1	1 1 1 1 1 1 1

- 7** Click the [Start] icon on the [Batch] assistant bar.  
The data is re-quantified using the corrected calibration curve.




### 5.3.3 Checking and Correcting Quantitation Results

Check the quantitation results for the samples with unknown concentrations.

- 1 Click the [Quantitative] icon on the [Postrun] assistant bar.



- 2 Double-click the data file to be checked from Data Explorer.  
The data file being checked opens.




- 3 Click the [Results] tab in the [Compound Table View].  
The quantitation results are displayed.

The image shows the 'Results' tab in the 'Compound Table View'. The table displays quantitation results for various compounds. The 'Conc.' column is mostly empty, indicating unknown concentrations.

ID#	Name	Type	ISTD G	m/z	Ret. Time	Ret. Index	Unit	Ref. Ions	Conc.
1	Naphthalene	ISTD	1	136.00	4.865	0	mg/L	108.00-13	
2	Acenaphthylene	Target	1	152.00	6.990	0	mg/L	151.00-15	0
3	Acenaphthen	ISTD	2	164.00	7.240	0	mg/L	162.00-16	
4	Fluorene	Target	2	166.00	8.200	0	mg/L	165.00-16	0
5	Phenanthrene	ISTD	3	188.00	10.250	0	mg/L	189.00-18	
6	Phenanthrene	Target	3	178.00	10.305	0	mg/L	176.00-17	0
7	Anthracene	Target	3	178.00	10.410	0	mg/L	176.00-17	0
8	Pyrene	Target	4	202.00	14.215	0	mg/L	200.00-20	0
9	Benzo[a]anthri	Target	4	228.00	17.790	0	mg/L	226.00-22	0
10	Chrysene-D12	ISTD	4	240.00	17.830	0	mg/L	236.00-24	
11	Chrysene	Target	4	228.00	17.895	0	mg/L	226.00-22	0
12	Benzo[b]fluor	Target	4	252.00	20.295	0	mg/L	250.00-25	0
13	Benzo[k]fluor	Target	5	252.00	20.340	0	mg/L	250.00-25	0
14	Benzo[a]pyre	Target	5	252.00	20.920	0	mg/L	250.00-25	0
15	Pyrene-D12	ISTD	5	264.00	21.035	0	mg/L	260.00-26	
16	Indeno[1,2,3-	Target	5	276.00	22.880	0	mg/L	274.00-27	0
17	Dibenz[a,h]peri	Target	5	276.00	22.915	0	mg/L	273.00-13	0
18	Benzo[ghi]peri	Target	5	276.00	23.330	0	mg/L	274.00-27	0


## 4 Display the standard spectra sub-window and reference data sub-window in the [Quantitative View] area.

If necessary, see ["Displaying Standard Spectra" P.72](#) , ["Displaying Reference Data" P.72](#) to display information about identified compounds.



## 5 Click on a compound name in the compound table and check the chromatogram in the [Quantitative View].

Check the results while viewing the peak identification/detection marks and baseline in the chromatogram.



### Reference

If necessary, perform identification or peak integration with reference to ["Manual Identification and Manual Peak Integration" P.63](#).


### NOTE

The same process can be accomplished more easily by performing the following operations on the chromatogram.

Process	Operation	Explanation
Manual Identification	[Shift] + [Ctrl] + right-click	Identifies integrated peaks.
Manual Peak Integration	[Shift] + right-click-drag	Connects start point and end point as baseline.
Manual Peak Integration	[Ctrl] + right-click-drag	Connects points with horizontal baseline.
Delete Identification Results	[Shift] + [Ctrl] + right-double-click	VOIDS identification and removes quantitative calculation results.

## 6 After checking the results, click (Save) on the toolbar.

The data file is saved.



### NOTE

When peaks are integrated for quantitation, concentrations calculated from the calibration curve are displayed.


However, if quantitation fails, character strings other than concentration values are displayed according to the cause of failure, as described below.

Displayed Character String	Explanation
No peak is detected.	Quantitative peak integration resulted in no peaks detected.
No peak is found in Window/Band range.	No peaks were detected within the retention time range specified for identification.
Ratio of reference ion does not match.	Peak is not identified due to the difference between specified and measured reference ion ratio values exceeding the allowable range.
Under the minimum similarity index.	Peak is not identified due to the measured similarity being less than the specified similarity setting, when mass pattern matching is specified in identification parameters.
No peak is identified.	Automatic identification results were manually deleted.

## ■ Displaying Standard Spectra

Data can be analyzed more easily by comparing the displayed spectrum with a standard spectrum.

- 1 Click [Spectrum View] on the [View] menu, then select [Display Standard Spectrum].  
The standard spectrum is displayed.



### NOTE

The standard spectrum is a mass spectrum of a standard sample registered when the compound table was created.

The standard spectrum can be hidden by repeating step 1 above.


When the measured spectrum is enlarged by dragging, the standard spectrum is enlarged correspondingly.

## ■ Displaying Reference Data

Compounds can be identified from the shape of chromatograms, retention times, and other information obtained by referencing measurement data of standard samples or spiked samples.


- 1 Select [Open Reference Data File] on the [File] menu to open the data file being referenced.

The reference data is displayed.



## 2 Select [Close Reference Data File] on the [File] menu to specify the reference data file to close.

The reference data file closes.



### NOTE

Up to three reference data files can be displayed.  
Reference data peaks cannot be integrated.


## 5.4 Printing Quantitative Analysis Reports

### 5.4.1 Creating and Outputting Quantitative Analysis Reports


#### 1 Click the [Report] icon on the [Quantitative] assistant bar. The [Data Report] window opens.



#### 2 Click [New Format File] on the [File] menu. The [File New] window opens.



### 3 Select [Use Template] and select the format [Quantitative Analysis Report].




### 4 Click [OK].

The [Quantitative Analysis Report] format opens.


### 5 Click the [Preview] icon on the [Data Report] assistant bar.

The print preview window opens.






# 6 After checking the report content, click [Print] to print the report.



5

# 7 Click (Save) on the toolbar. The report is saved as a data file.



# 6


## Shutting Down GC/MS

### 6.1 Vacuum System Shutdown


- 1 Click the [Vacuum Control] icon on the [Real Time] assistant bar.  
The [Vacuum Control] window opens.



- 2 Click [Auto Shutdown].  
The vacuum system shuts down.



- 3 When [Completed] is displayed, click [Close].



## 6.2 Turning OFF the Power

Turn OFF the power by performing the procedure for turning ON the power in reverse.







If accessory/peripheral equipment is connected, switch OFF the accessory/peripheral equipment power last.

### Reference

See "[2.1 Turning ON the Power](#)" P.3 for details on how to turn ON the power.

- 1** Quit the [GCMS Real Time Analysis] program and all other programs that are running.
- 2** Turn OFF the power to the PC, printer, and display.
- 3** Turn OFF the power to the MS unit.
- 4** Turn OFF the power to the GC unit.

GCMSsolution uses the file formats described below.

File type	Icon	Extension	File contents
Data file		.qgd	In addition to the raw data acquired (e.g., chromatograms and spectra), the following information is saved. <ul style="list-style-type: none"> <li>• Calculation results such as area values and concentrations</li> <li>• Status information such as the oven temperature and error status at the time data is acquired</li> <li>• Contents of method files used in analysis (including configuration settings used for analysis)</li> <li>• Contents of report format file (when reports are output)</li> <li>• Contents of batch files (when batch processing is performed)</li> <li>• Contents of tuning file used in analysis</li> </ul>
Method file		.qgm	Analysis conditions, peak integration parameters, compound tables, etc. are saved. Because the configuration settings are saved when the method is edited, the configuration settings are checked when the method file is loaded to ensure that they agree with the current settings. Created calibration curves are also saved in the method file.
Report format file		.qgr	The report format information used to output a report, such as layout information and detailed settings, is saved. Once a report format file has been created, it can be used repeatedly to output reports of the same format.
Batch file		.qgb	Batch tables used to perform automatic sequential processing are saved. The same files can be used in both the [GCMS Real Time Analysis] program and the [GCMS Postrun Analysis] program.
Tuning file		.qgt	The conditions used to perform instrument adjustment (tuning) and the tuning results are saved.
Library file		.lib	These files are used to register the compound information and spectral data used to perform similarity searches. The libraries consist of public libraries (e.g., NIST and Wiley) and private libraries.

# Appendix **B**

## Viewing Help

### B.1 Viewing Help

If you do not know how to perform a procedure, refer to Help using one of the procedures described below.


#### B.1.1 Displaying Help from the Assistant Bar

- 1 Click the [Guide] icon on the assistant bar in a window to display explanations for that window.



#### B.1.2 Displaying Help from the Menu Bar

- 1 Click [Contents] on the [Help] menu displayed in the menu bar of a window to display the [GCMS Help window].

**B**

#### Searching from the [Contents] Tab

- 1 Double-click the applicable topic.

### Searching from the [Index] Tab

- 1 Type the applicable word.  
The topic that matches the word is displayed at the top.
- 2 Select the applicable topic and click [Display].  
Details on the selected topic are displayed.



#### NOTE

If there are multiple instances of the word, the [Topics Found] window opens. Select the applicable topic from the list displayed and click [Display].

### Searching from the [Search] Tab

- 1 Type the applicable word and click [Search].  
The search results are displayed.
- 2 Select the applicable topic and click [Display].  
Details on the selected topic are displayed.

## B.1.3 Displaying Help with the F1 Key

**1**


**Press the [F1] key on the keyboard.**



Help for the open window is displayed.

# Appendix C

## Using Data Explorer

Using the data explorer function makes it easy to load files.  
For convenience, leave Data Explorer displayed at all times.




No.	Explanation
1	Click  (Data Explorer) on the toolbar to open and close the [Data Explorer] window.
2	Click  (Project (Folder) Selection) to open the [Project (Folder) Selection] window and then click the applicable folder.
3	Click [Close] to close the [Project (Folder) Selection] window.
4	Click a tab to display the file of the corresponding type.
5	Click a file name (6) to display the corresponding file information.
6	Double-click the appropriate file to load the data.



**NOTE**


[Create New Project (Folder)] can be used to copy files from the folder currently open in Data Explorer and to create new folders.



This creates a new folder at the same directory level as currently open in Data Explorer.

Using this feature is especially useful in the following situation.

Example:






# Appendix **D** Maintenance

## D.1 Maintenance

Replace or clean the consumable items and maintenance parts as necessary, referring to the [MS Navigator] window using the procedure described below.

- 1** Double-click the  (GCMS Real Time Analysis) icon.  
The [GCMS Real Time Analysis] program starts.

- 2** Select [Maintenance] on the [Help] menu.  
The [MS Navigator] window opens.




- 3** Click on the instrument for which maintenance will be performed.




# 4

Read the precautionary information carefully and then click the applicable item under the maintenance menu.



## 5 Perform maintenance by following the instructions displayed on the screen.

Click an image to enlarge it. Click [Back] in the enlarged window to return to the original window.



To perform another maintenance item, click [Back] and repeat the procedure from step 3.

After completing maintenance, close the [MS Navigator] window.

After performing maintenance, reset the usage frequencies and usage times using the procedure described in ["Appendix D.3 Reset Procedure for Usage Frequencies and Usage Times" Ver.2.5](#) P.87.

## D.2 Easy sTop (Applicable to QP2010 Ultra and QP2010 SE models)

Using Easy sTop allows replacing septa and glass inserts without stopping the vacuum system. Therefore, it significantly reduces the time required for stabilizing the system after replacement and eliminates the need for autotuning.

To protect columns, Easy sTop keeps the temperature of the sample injection unit, column oven, and interface at 70 °C or below. Consequently, it can take about 30 minutes until glass inserts and septa can be replaced.


### 1 Double-click one of the icons for consumables in the instrument monitor.

The [Consumable] tab page opens in the [Monitor Settings] window.




## 2 Click [Easy sTop].

The [Easy sTop] window opens and the injection unit, column oven, and interface temperatures decrease. When each temperature reaches 70 °C or lower, the "Ready" status is displayed in the [Easy sTop] window.



## 3 Click [Replace], then replace septa or glass inserts in the sample injection unit.

For replacement procedures, refer to the septum replacement procedure or insert replacement procedure in the [MS Navigator] window.




 **NOTE**

Clicking [Replace] stops the supply of carrier gas.

If left in that state for extended periods, it could reduce column performance. Therefore, replace septa and inserts as quickly as possible.

**4** After replacement, click [Complete] in the [Easy sTop] window.

If there is no air leaking in, the sample injection unit, column oven, and interface temperatures return to their previous temperatures before Easy sTop started.

**5** Reset the usage counter for the septum and glass insert.

For instructions on how to reset usage counters, see the procedure on page 88, starting with step 3.

## D.3 Reset Procedure for Usage Frequencies and Usage Times Ver.2.5

The GCMSsolution's instrument monitor keeps track of the usage frequencies and usage times of the consumable items.

After replacing a consumable item or cleaning the ion source, reset the usage frequencies and usage times using the procedure described below.


**1** Double-click one of the icons for consumables in the instrument monitor.

The [Consumable] tab page opens in the [Monitor Settings] window.




D

- 2 Click **[Reset Consumables]**.  
The **[Reset Consumables]** window opens.



- 3 Click **[Clear]** for the consumable items subjected to maintenance or input "0" in the usage time cells as appropriate, and click **[OK]**.  
The previous window returns.



- 4 Click **[OK]** in the **[Monitor Settings]** window.  
The window closes.

## D.4 Changing Replacement Guidelines for septa and Glass Inserts


For septa, replacement frequency varies depending on the syringe needle diameter.

Glass insert replacement frequency varies depending on the sample. Set replacement guidelines based on the sample.

- 1 Click the [System Configuration] icon on the [Real Time] assistant bar.  
The [System Configuration] window opens.




- 2 Double-click [SPL1] under [Modules Used for Analysis].  
The [Modules of Analytical Line #1] window opens.




### 3 Click [Injection Port Maintenance].

The [Injection Port Maintenance (SPL1)] window opens.



### 4 Input [Septum Used Counts] and [Insert Used Counts] settings.

To restore default settings, click [Default].



### 5 Click [OK].


The [Modules of Analytical Line #1] window returns.

### 6 Click [OK].

The [System Configuration] window returns.



# 7 Click [Set].



The replacement guidelines for septa and glass inserts are changed.



# Single Analysis (Manual Injection)


Use the procedure described below when analyzing samples one-by-one using the autosampler or when performing analysis using manual injection.

- 1** Start the [GCMS Real Time Analysis] program, then click the [Data Acquisition] icon on the [Real Time] assistant bar.

The [Acquisition] window opens.




- 2** Double-click the method file to be used in Data Explorer.



- 3** Click the [Sample Login] icon on the [Acquisition] assistant bar.



The [Sample Login] window opens.



- 1 Enter [Sample Name] and [Data File].
- 2 When using an autosampler, input [Vial #] in which the sample is set and [Injection Volume].
- 3 Click [OK].


#### 4 When using an autosampler, set syringe rinse solvent and samples in the specified positions.

#### 5 Click the [Download] icon on the [Acquisition] assistant bar.

The method file settings are transferred to the instrument. When preparation for GC and MS has been completed, the [Start] icon turns green, indicating that it can be selected. If using autosampler model AOC-20i, the analysis starts automatically.



#### 6 Inject the sample and press the [START] button on the keyboard at the GC unit.

If using accessory/peripheral equipment, start such equipment first, then click the  (Start) icon.

#### NOTE

To abort analysis before completion, click the  (Stop) icon on the [Acquisition] assistant bar.

## Index Searches

It is possible to search for information related to the target compounds (e.g., spectra and information on the structure) in the library.


### 1 Click the [Library Editor] icon on the [Postrun] assistant bar.

The [Library Editor] window opens.



### 2 Click [Open Library] on the [File] menu to open the library to be used.

The library opens.



### 3 Click the second row cell in the [Index] column to select an item.

Parameters		Parameter	Upper/Lower	# of Hit
1	Index	1-191436	No need to set	191436
2	No Setting		No need to set	0
3	Serial Number		No need to set	0
4	Mol Wt		No need to set	0
5	Formula		No need to set	0
6	Cmpd Name		No need to set	0
7	Base Peak		No need to set	0
8	Ret. Index		No need to set	0
9	Class Flags		No need to set	0
10	CAS Number		No need to set	0
11	No Setting		No need to set	0

### 4 Enter information for the index item in the [Parameter] column for the row where the index item was selected.

Parameters		Parameter	Upper/Lower	# of Hit
1	Serial Number	1-191436	No need to set	191436
2	Cmpd Name	hexa	<input type="checkbox"/> Match Case	0
3	No Setting		No need to set	0
4	No Setting		No need to set	0
5	No Setting		No need to set	0
6	No Setting		No need to set	0

## 5 Click [Start] on the [Index Search] menu.

The results are displayed.

Add index items until the target compound is found.

n Analysis (Admin) - [Library Editor - NIST08.LIB (191,436 Spectrum)]

List Target Info. Compound Info. Index Search Tools Window Help

Start  
Export Search Results

Index	Parameter
1 Serial Number	1-191436
2 Cmpd Name	hexa
3 Formula	C6
4 Cmpd Name	benzene

## 6 Confirm the applicable information (e.g., spectrum or structure).

Parameters

Index	Parameter	Upper/Lower	# of Hit
1 Serial Number	1-191436	No need to set	131436
2 Cmpd Name	hexa	<input type="checkbox"/> Match Case	11845
3 Formula	C6	<input checked="" type="checkbox"/> Match Case	146
4 Cmpd Name	benzene	<input type="checkbox"/> Match Case	13
5 Mol Wt	282-282	No need to set	1
6 No Setting		No need to set	0

Hit#	Cmpd Name	Mol Wt	Formula
1	Benzene, hexachloro- Perchlorobenzene	282	C6Cl6

1: 282: Benzene, hexachloro- Perchlorobenzene Amatin Anticane Bunt-cure Bunk-no more Co-op Hexa HCB Juin's Carbon chloride No Bunt No Bunt Liquid No Bunt 40 No Bunt 80 Pentachlorophenyl chlorid

Inten. (x10,000)

Base Peak: 284/10,000

C1=CC=C(C=C1)Cl

Mass spectrum plot showing intensity vs m/z. Major peaks at 107, 142, 144, 177, 214, 249, 284.

CAS#: 118-74-1 Mol Wt 282 Serial#: 98409  
Cmpd Name: Benzene, hexachloro- Perchlorobenzene Amatin Anticane Bunt-cure Bunk-no more Co-op Hexa HCB Juin's Carbon chloride No Bunt No Bunt Liquid No Bunt 40 No Bunt 80 Pentachlorophenyl c  
Cmpd Form: C6Cl6 Class Flag No Class Flags  
Description:  
Ret Index: 1760

# Displaying Chromatograms

Displaying the appropriate mass chromatogram while analyzing data for qualitative analysis makes analysis easier.

If the  $m/z$  value for which the chromatogram is to be displayed is not known, search for it beforehand using the procedure described in ["Appendix F Index Searches" P.94.](#)


## G.1 Displaying Chromatograms from Fragment Tables

- 1 Click the [Fragment Table] icon on the [Qualitative] assistant bar.  
The [MC Fragment Table] window opens.



- 2 Enter the applicable values in the [ $m/z$ ] and [Factor] columns, select the corresponding cells in the [Disp.] column, and click [OK].

A mass chromatogram is displayed in the MC window.




The display can be changed as shown below by enabling/disabling [Base Shift] in the table.

- With Base Shift




- Without Base Shift




## G.2 Displaying Chromatograms from Mass Spectra

- 1 In the mass spectrum, specify and enlarge the range containing the desired peaks by dragging the mouse.



- 2 Move the mouse pointer to the spectral peak to be displayed and double-click. A mass chromatogram is displayed in the MC window, enlarged by an automatically set enlargement rate.



### NOTE

- To hide the mass chromatogram, deselect the applicable cell in the [Disp.] column in the [MC Fragment Table] window.
- To undo enlarging, right-click on the mass spectrum and select [Undo Zoom] on the menu that appears.



# Appendix H

## Editing and Adding Batch Files During Sequential Analysis

During sequential analyses, the batch file being executed can be edited while the analysis is still in progress. In addition, a separate batch file can be added for processing after the current batch file is finished, by specifying a batch queue.

### H.1 Editing Batch Files

- 1 Click on the batch table, then after the window switches, click the [Pause/Restart] icon on the [Batch] assistant bar.

The [Batch Table] window opens, allowing unexecuted rows to be edited.

The screenshot displays the GCMS Real Time Analysis software interface. At the top, a table lists batch rows with columns for View#, Sample Name, Sample ID, Sample Type, Analysis Type, Method File, Data File, Level#, Inj. Volume, ISTD Amt, Report Output, Report File, and Tuning File. Below the table, a 'Batch' window is open, showing a file explorer view of the batch files. A 'Line1' window is also open, displaying a chromatogram and analysis parameters. A text box overlay on the chromatogram reads: 'When analyses are in progress, both the [Batch Table] window and [Acquisition] window are displayed. To switch to the [Batch Table] window, click on the table.'

View#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Level#	Inj. Volume	ISTD Amt	Report Output	Report File	Tuning File
1	Methylene Chloride	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank001.qgd	1	1	1	[Level1 Con]		
2	Standard Sample	STD-0001	1:Standard (I)	IT QT	PAH_SIM.qgm	20091115_STD5.qgd	1	1	1	[Level1 Con]		
3	Standard Sample	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD6.qgd	2					
4	Standard Sample	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD7.qgd	3					
5	Standard Sample	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD8.qgd	4					
6	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample3.qgd	1					
7	Unknown Sample	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample4.qgd	1					




#### NOTE

Analysis of rows currently being analyzed will continue to be executed.


## 2 Edit the batch table.

Right-click on the row to be edited, then select [Add Row], [Delete Row], or other action on the menu that appears.

The vial number, data file name, or other information can be changed as well.



## 3 Click (Save) on the toolbar.



## 4 Click the [Pause/Restart] icon on the [Batch] assistant bar. The analysis restarts.



### NOTE

Some accessory/peripheral equipment may prevent using this function.

## H.2 Adding Batch Files (Batch Queue)

### H.2.1 Creating Batch Files to Add

**1** Double-click the  (GCMS Analysis Editor) icon.

The [GCMS Analysis Editor] program starts.

**2** Click the [Batch Processing] icon on the [Real Time] assistant bar.

The [Batch Table] window opens.




**3** Create the batch file to be added.

The screenshot shows the GCMS Analysis Editor (Admin) - [Batch Table - Untitled] window. The main area displays a table with the following data:

Visit#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Level#	Inj. Volume	ISTD
1	Unknown Sample	UNK-0001	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample3.ogd	1	1	(Level1
2	Unknown Sample	UNK-0002	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample4.ogd	1	1	(Level1
3	Unknown Sample	UNK-0003	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample5.ogd	1	1	(Level1
4	Unknown Sample	UNK-0004	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample5.ogd	1	1	(Level1
5	Unknown Sample	UNK-0005	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample7.ogd	1	1	(Level1
6	Unknown Sample	UNK-0006	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample8.ogd	1	1	(Level1
7	Unknown Sample	UNK-0007	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample9.ogd	1	1	(Level1
8	Unknown Sample	UNK-0008	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample10.ogd	1	1	(Level1
9	Unknown Sample	UNK-0009	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample11.ogd	1	1	(Level1
10	Unknown Sample	UNK-0010	Q:Unknown	IT QT	PAH_SIM.agm	20091115_sample12.ogd	1	1	(Level1

## 4 Name and save the batch file.



## 5 Quit the [GCMS Analysis Editor] program.



### NOTE

- The analysis will not start if the same data file name is used more than once or the specified method file does not exist.
- The batch queue is not activated until the [GCMS Analysis Editor] program is closed.

## H.2.2 Adding Batch Files

### 1 Start the [GCMS Real Time Analysis] program.

During analysis, both the [Acquisition] and [Batch Table] windows are displayed simultaneously.

**[Batch Table] Window**

Line#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Level#	Inj. Volume	ISTD Amt.	Report Output	Report File	Tuning File
1	Methylene Chloride	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank001.qgd	1	1	Level1 Con	Print		
2	Standard Sample	STD-0001	1:Standard (I)	IT QT	PAH_SIM.qgm	20091115_STD5.qgd	1	1	Level1 Con	Print		
3	Standard Sample	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD6.qgd	2	1	Level1 Con	Print		
4	Standard Sample	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD7.qgd	3	1	Level1 Con	Print		
5	Standard Sample	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD8.qgd	4	1	Level1 Con	Print		
6	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample3.qgd	1	1	Level1 Con	Print		
7	Unknown Sample	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample4.qgd	1	1	Level1 Con	Print		

**[Acquisition] Window**

Sample Name: Methylene Chloride  
 Sample ID: UNK-0001  
 Data Description: Line1-MS Ret. Time: 9.867 / 25.020 min Scan #: 1761  
 Max Intensity: 12,989

Chromatogram showing intensity vs. time (min) with peaks at 1.172, 1.172, and 1.172 min.

**Data Explorer - Method**

Project in: ..\20091115\_operation guide

File Name: PAH\_Scan (11/15/2009 4:30 PM), PAH\_SIM (11/13/2009 8:05 AM)

Method Parameters:

- Inj. Port: SPL1 Inj. Heat Port: INJ1
- Column Oven Temp.: 45.0 °C
- Injection Temp.: 250.0 °C
- Injection Mode: Splitless
- Sampling Time: 1.00 min
- Carrier Gas: He Pim. Press.: 900-900
- Flow Control Mode: Linear Velocity
- Pressure: 67.7 kPa
- Total Flow: 50.0 mL/min
- Column Flow: 1.22 mL/min
- Linear Velocity: 40.0 cm/sec
- Purge Flow: 3.0 mL/min


Program: Column Oven Temperature

Rate	Final Temperature	Hold Time
0	45.0	1.00
1	130.0	0.00
2	180.0	0.00
3	240.0	0.00

Total Program Time: 26.29 min


## 2 Click the [Batch Table] window.

The content of the toolbar, menu bar, and assistant bar changes.




## 3 Select [Batch Queue] on the [Batch] menu.

The [Batch Queue] window opens.




**4** Click [Add] to open the batch file to be added.



**5** If multiple batch files were added, change their order by clicking to select the desired batch file, then clicking [Move Up] or [Move Down].

Files are executed consecutively from the top.



**6** When finished editing, click [OK].

# Reducing the Carrier Gas Flow Rate After Sequential Analysis and Using the Ecology Mode

Reducing the carrier gas flow rate after analysis is finished is recommended to reduce carrier gas consumption.

## I.1 Reducing the Carrier Gas Flow Rate After Sequential Analysis

For models other than QP2010 Ultra and SE, perform the following operations.

### I.1.1 Creating a Method File That Reduces the Carrier Gas Flow Rate

As an example, the following describes how to create a method file that reduces the total flow rate to 20 mL/min.



- 1 Start the [GCMS Real Time Analysis] program, then in Data Explorer, double-click the method file to be used for sequential analysis.

The screenshot displays the GCMS Real Time Analysis software interface. The main window shows a method configuration for 'PAH\_SIM'. The 'Column Oven Temperature' program is selected, and its parameters are shown in the 'Program' table below.

Rate	Final Temperature	Hold Time
0	45.0	1.00
1	130.0	0.00
2	180.0	0.00
3	240.0	0.00

The 'Total Program Time' is 26.29 min. The 'Column' details are: Name: Rtx-SMS, Length: 30.0 m, Thickness: 0.25 μm, Diameter: 0.25 mm. The 'Detail of Injection Port...' section shows 'High Press. Injection' and 'Carrier Gas Saver' are selected. The 'Ready Check...' section shows 'GC Program...' is selected.

## 2 Change [Total Flow] to 20 mL/min, then name and save the method file.




### I.1.2 Creating Batch Files

- 1 Click the [Batch Processing] icon on the [Real Time] assistant bar. The [Batch Table] window opens.




- 2 In Data Explorer, double-click the batch file to be used for sequential analysis.






- 3** Right-click on the batch table and select [Table Style] on the menu that appears. The [Table Style] window opens.



- 4** Click [Run Mode] in the [Hide Items] list, then click [Add>>] and [OK]. A [Run Mode] column is added to the end of the batch schedule.



## 5 Edit the batch file.

Add a row at the end and select a method file created in "Appendix I.1.2 Creating Batch Files" P.106. Vial number, level number, and injection volume settings do not need to be changed from their default values. Enter a data file name that is not the same as any other row.


Folder: C:\GCMSsolution\20091115_operation guide							
	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File
1	1	Methylene Chloride	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd
2	2	Standard Sample	STD-0001	1:Standard:(f)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd
3	3	Standard Sample	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd
4	4	Standard Sample	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd
5	5	Standard Sample	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd
6	1	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample1.qgd
7	6	Unknown Sample	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample2.qgd
8	1			0:Unknown	IT QT	PAH_SIM_low.qgm	20091115_low.qgd

## 6 Click the [Run Mode] cell for the row that specifies the method file that reduces the flow rate, then click the arrow button that appears.

The [Run Mode] window opens.

Folder: C:\GCMSsolution\20091115_operation guide								
	Level#	Inj. Volume	ISTD Amt.	Report Output	Report File	Tuning File	Data Description	Run Mode
1	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
2	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
3	2	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
4	3	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
5	4	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
6	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
7	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
8	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP
9	1	1	(Level1 Con	<input type="checkbox"/> Print				DL AQ DP

## 7 Configure [Run Mode] settings as shown below, then click [OK].



## 8 Name and save the batch file, then click the [Start] icon on the [Batch] assistant bar.

This results in loading the method file that reduces flow rate after the analysis for the seventh row is finished, which ends the sequential analysis with the carrier gas flow rate at 20 mL/min.



## I.2 Ecology Mode (This feature applies to QP2010 Ultra and QP2010 SE models.)


Using the ecology mode reduces power consumption and carrier gas consumption during standby for analysis.

### I.2.1 Setting the Mode Manually

- 1 Click the [Ecology Mode] icon in the instrument monitor.  
A message window opens.



- 2 Click [Yes].  
The [Ecology Mode] window opens and the mode switches to the ecology mode.  
After switching to the ecology mode, the column oven temperature and the total carrier gas flow rate decrease.




#### NOTE

The [Ecology Mode] window is displayed when in the ecology mode. Cancel the ecology mode before using [GCMS Real Time Analysis] to perform operations in other windows.

To cancel the ecology mode, click [Cancel] in the [Ecology Mode] window.

When the ecology mode is canceled, settings before switching to the ecology mode are restored.



## 1.2.2 Setting the Mode Using Batch Processing

This allows switching the instrument to the ecology mode after the entire sequential analysis is finished.


- 1 Click the **[Batch Processing]** icon on the **[Real Time]** assistant bar.  
The **[Batch Table]** window opens.



- 2 Create and save a batch file.

Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Len
1	Methylene Chloride	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_blank.qgd	
2	Standard Sample	STD-0001	1:Standard (I)	IT QT	PAH_SIM.qgm	20091115_STD1.qgd	
3	Standard Sample	STD-0002	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD2.qgd	
4	Standard Sample	STD-0003	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD3.qgd	
5	Standard Sample	STD-0004	1:Standard	IT QT	PAH_SIM.qgm	20091115_STD4.qgd	
6	Blank	UNK-0001	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample1.qgd	
7	Unknown Sample	UNK-0002	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample2.qgd	
8	Unknown Sample	UNK-0003	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample3.qgd	
9	Unknown Sample	UNK-0004	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample4.qgd	
10	Unknown Sample	UNK-0005	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample5.qgd	
11	Unknown Sample	UNK-0006	0:Unknown	IT QT	PAH_SIM.qgm	20091115_sample6.qgd	

**3** Select [Enters Ecology Mode when ending Realtime Batch] on the [Batch] menu.




**4** Click the [Start] icon on the [Batch] assistant bar.



**5** When the ecology mode confirmation message appears, click [Yes].

The mode switches to the ecology mode after the sequential analysis is completely finished, including the batch queue.




**NOTE**

The setting can be canceled by repeating step 3, but leave the setting as it is.

The [Ecology Mode] window is displayed when in the ecology mode. Cancel the ecology mode before using [GCMS Real Time Analysis] to perform operations in other windows.

To cancel the ecology mode, click [Cancel] in the [Ecology Mode] window.

When the ecology mode is canceled, settings before switching to the ecology mode are restored.




# Printing Reports

Reports can be output from GCMSsolution using the two methods described below.


- Image printing : The image in the displayed window is automatically converted to a report.
- Report creation : A report format is set and output manually.

## J.1 Printing Images (Printing Spectra and Chromatograms Displayed in Windows)


- 1 Call up the applicable data in the [Data Analysis] window in the qualitative or quantitative processing modes of the [GCMS Postrun Analysis] program.
- 2 Display the chromatogram and mass spectrum in the window in the way desired for the report.



- 3 Point to [Print Image] on the [File] menu and select [Edit Format].  
The [Report] window opens.



- 4 Adjust the size as necessary.



- 5** After editing, click the [Print] icon on the [Report] assistant bar.  
The report is output.




- 6** After outputting the report, close the [Report] window.

## J.2 Creating Reports

With report creation, reports are output after setting report formats or using previously created templates.


Process and save the results to be output (such as spectral information) in advance.

- 1** Open the applicable data in the [GCMS Postrun Analysis] - [Data Analysis] window.  
The same report is output for both the qualitative and quantitative windows.

- 2** Click the  (Report) icon on the [Qualitative] or [Quantitative] assistant bar.  
The [Data Report] window opens.


### J.2.1 Using Templates

- 1** Select [New Format File] on the [File] menu.





- 2** Select [Use Template], select the applicable template, and click [OK].




**NOTE**

If this selection window is not displayed, select [Option] on the [Tool] menu to display the [Setting Options] window and, on the [File New] tab, select [Prompt on File New] for the report format file.












## J.2.2 Using Previously Created Report Files

- 1** In Data Explorer, double-click the report file to be used.




## J.2.3 Manually Setting Report Content


- 1 Click the buttons on the toolbar for the information to be printed or select the desired items on the [Item] menu.

Icon	Name	Explanation
	Sample information	Select to print sample information.
	Method	Select to print methods.
	Peak table	Select to print the peak tables in qualitative tables.
	Chromatogram	Select to print the chromatograms (TIC, MIC, and MC).
	Spectrum graph	Select to print the mass spectra registered in spectrum processing tables.
	Mass table	Select to print the mass tables for the spectra registered in spectrum processing tables.
	Quantitative graph	Select to print the chromatograms and quantitative values obtained in quantitative results.
	Quantitative table	Select to print the tables obtained in quantitative results.
	Calibration curve	Select to print calibration curves.
	Tuning	Select to print the tuning results obtained when data acquisition is executed.
	Library search	Select to print the library search results obtained for the mass spectra registered in spectrum tables. <ul style="list-style-type: none"> <li>• Searches must be performed in the spectrum tables.</li> </ul>

- 2 Drag the mouse in the layout view to specify the print range.**  
The properties window for the item being laid out opens.



- 3 Set [Properties] and click [OK].**




**Reference**

Refer to Help for details on property settings.


**NOTE**

To display a properties window again, double-click on the corresponding item.

- 4** Click the [Preview] icon on the [Data Report] assistant bar and check the contents of the report being output.



- 5** After the checking the report content, click [Print] to output the report.



- 6** Select [Save Format File As] on the [File] menu to name and save the report file. This allows loading the report format in the future to create reports easily.

